

aks-saniflow III



Alternating pressure system
of the middle standard class
for pressure sore
up to grade III
with static function



- stageless pressure adjustment
- save plug-in connection
- optical pressure gage with 2 LED
- acoustic alarm when pressure is low
- no risk of confusion
with the mattress labelling



Stand: 1/2007

aks

aktuelle krankpflege systeme

Antwepener Straße 6
53842 Troisdorf
Fon: 02241/94 74-0
Fax: 02241/94 74-88
E-mail: aks@aks.de
Web: www.aks.de



aks-saniflow III

Alternating pressure system of the middle standard class for pressure sore up to grade III with static function

Pressure pump

- nice and functional design
- cycle time: 12 minutes
- pressure adjustment according to the patient's weight (picture 1)
- switch to deactivate the alarm (picture 2)
- switch for static function (picture 2)
- stable and retractable metal handles
- easy filter change
- very elastic foam strap and soft feet on the housing to avoid any noise



Mattress

- overlay system
- alternating 2-chamber system for alternating pressure relief
- 17 cells
- all cells are independently interchangeable
- cell height: 13 cm
- breathable and waterproof NYLON/PVC complete cover closed with zipper and fixed on the mattress
- diagonal rubber bands over the edges for fixing on the bed mattress
- quick-air-removing device (picture 3)



Service

- competent service check of the pump for the re-installation at fixed price (incl. filter change, without repairs and without spare parts, with electrical testing after DIN VDE 0751-1)
- repair of pump at a fixed price incl. spare parts
- hygiene and transport box for safe transport
- training at aks in Germany or when required on-site

Rights reserved to alter specifications

Technical data

Pump

Electrical power supply	230 Volt, 50 Hz, 12 W
Weight	2,4 kg
Dimensions	W 28 x H 10 (incl. feet) x L 20,5 cm
Warranty	2 years

Mattress

Dimensions	90 x 200 cm
Height of cells	13 cm
Quantity of cells	17
Patient's weight	up to 120 kg

Warranty 1 year



aks-Hygiene- and transportation box for closed hygiene cycle.



Ref	aks-saniflow III	S21233
	Replacement pump	S21234
	Replacement mattress	S21235